








Tube #	Salt	Tube #	Polymer	Tube #	pH $\diamond$	
1.	0.2 M Sodium fluoride	1.	20% w/v Polyethylene glycol 3,350	1.	7.3	F <sup>-</sup>
2.	0.2 M Potassium fluoride	2.	20% w/v Polyethylene glycol 3,350	2.	7.3	Cl <sup>-</sup>
3.	0.2 M Ammonium fluoride	3.	20% w/v Polyethylene glycol 3,350	3.	6.2	I <sup>-</sup>
4.	0.2 M Lithium chloride	4.	20% w/v Polyethylene glycol 3,350	4.	6.8	Fluoride Chloride Iodide
5.	0.2 M Magnesium chloride hexahydrate	5.	20% w/v Polyethylene glycol 3,350	5.	5.9	
6.	0.2 M Sodium chloride	6.	20% w/v Polyethylene glycol 3,350	6.	6.9	
7.	0.2 M Calcium chloride dihydrate	7.	20% w/v Polyethylene glycol 3,350	7.	5.1	
8.	0.2 M Potassium chloride	8.	20% w/v Polyethylene glycol 3,350	8.	7.0	
9.	0.2 M Ammonium chloride	9.	20% w/v Polyethylene glycol 3,350	9.	6.3	
10.	0.2 M Sodium iodide	10.	20% w/v Polyethylene glycol 3,350	10.	7.0	
11.	0.2 M Potassium iodide	11.	20% w/v Polyethylene glycol 3,350	11.	7.0	
12.	0.2 M Ammonium iodide	12.	20% w/v Polyethylene glycol 3,350	12.	6.2	
13.	0.2 M Sodium thiocyanate	13.	20% w/v Polyethylene glycol 3,350	13.	6.9	
14.	0.2 M Potassium thiocyanate	14.	20% w/v Polyethylene glycol 3,350	14.	7.0	
15.	0.2 M Lithium nitrate	15.	20% w/v Polyethylene glycol 3,350	15.	7.1	
16.	0.2 M Magnesium nitrate hexahydrate	16.	20% w/v Polyethylene glycol 3,350	16.	5.9	
17.	0.2 M Sodium nitrate	17.	20% w/v Polyethylene glycol 3,350	17.	6.8	
18.	0.2 M Potassium nitrate	18.	20% w/v Polyethylene glycol 3,350	18.	6.8	
19.	0.2 M Ammonium nitrate	19.	20% w/v Polyethylene glycol 3,350	19.	6.2	
20.	0.2 M Magnesium formate dihydrate	20.	20% w/v Polyethylene glycol 3,350	20.	7.0	
21.	0.2 M Sodium formate	21.	20% w/v Polyethylene glycol 3,350	21.	7.2	
22.	0.2 M Potassium formate	22.	20% w/v Polyethylene glycol 3,350	22.	7.3	
23.	0.2 M Ammonium formate	23.	20% w/v Polyethylene glycol 3,350	23.	6.6	
24.	0.2 M Lithium acetate dihydrate	24.	20% w/v Polyethylene glycol 3,350	24.	7.9	
25.	0.2 M Magnesium acetate tetrahydrate	25.	20% w/v Polyethylene glycol 3,350	25.	7.9	
26.	0.2 M Zinc acetate dihydrate	26.	20% w/v Polyethylene glycol 3,350	26.	6.4	
27.	0.2 M Sodium acetate trihydrate	27.	20% w/v Polyethylene glycol 3,350	27.	8.0	
28.	0.2 M Calcium acetate hydrate	28.	20% w/v Polyethylene glycol 3,350	28.	7.5	
29.	0.2 M Potassium acetate	29.	20% w/v Polyethylene glycol 3,350	29.	8.1	
30.	0.2 M Ammonium acetate	30.	20% w/v Polyethylene glycol 3,350	30.	7.1	
31.	0.2 M Lithium sulfate monohydrate	31.	20% w/v Polyethylene glycol 3,350	31.	6.0	
32.	0.2 M Magnesium sulfate heptahydrate	32.	20% w/v Polyethylene glycol 3,350	32.	6.0	
33.	0.2 M Sodium sulfate decahydrate	33.	20% w/v Polyethylene glycol 3,350	33.	6.7	
34.	0.2 M Potassium sulfate	34.	20% w/v Polyethylene glycol 3,350	34.	6.8	
35.	0.2 M Ammonium sulfate	35.	20% w/v Polyethylene glycol 3,350	35.	6.0	
36.	0.2 M Sodium tartrate dibasic dihydrate	36.	20% w/v Polyethylene glycol 3,350	36.	7.3	
37.	0.2 M Potassium sodium tartrate tetrahydrate	37.	20% w/v Polyethylene glycol 3,350	37.	7.4	
38.	0.2 M Ammonium tartrate dibasic	38.	20% w/v Polyethylene glycol 3,350	38.	6.6	
39.	0.2 M Sodium phosphate monobasic monohydrate	39.	20% w/v Polyethylene glycol 3,350	39.	4.7	
40.	0.2 M Sodium phosphate dibasic dihydrate	40.	20% w/v Polyethylene glycol 3,350	40.	9.1	
41.	0.2 M Potassium phosphate monobasic	41.	20% w/v Polyethylene glycol 3,350	41.	4.8	
42.	0.2 M Potassium phosphate dibasic	42.	20% w/v Polyethylene glycol 3,350	42.	9.2	
43.	0.2 M Ammonium phosphate monobasic	43.	20% w/v Polyethylene glycol 3,350	43.	4.6	
44.	0.2 M Ammonium phosphate dibasic	44.	20% w/v Polyethylene glycol 3,350	44.	8.0	
45.	0.2 M Lithium citrate tribasic tetrahydrate	45.	20% w/v Polyethylene glycol 3,350	45.	8.4	
46.	0.2 M Sodium citrate tribasic dihydrate	46.	20% w/v Polyethylene glycol 3,350	46.	8.3	
47.	0.2 M Potassium citrate tribasic monohydrate	47.	20% w/v Polyethylene glycol 3,350	47.	8.3	
48.	0.2 M Ammonium citrate dibasic	48.	20% w/v Polyethylene glycol 3,350	48.	5.1	



$\diamond$  Measured pH at 25 ° C

PEG/Ion Screen contains forty-eight unique reagents. To determine the formulation of each reagent, simply read across the page.